


Trace Element Analysis of Industrial Wastewater and Sewage with TXRF



Bruker Nano GmbH
Berlin, Germany



Welcome



Today's Topics

- TXRF – how does it work?
- Wastewater analysis
- Water monitoring with bioindicators
- Comparison with Atomic Spectroscopy methods
- Interactive Q & A

Speakers

Dr. Hagen Stosnach
Applications Scientist TXRF
Berlin, Germany



Dr. Armin Gross
Global Product Manager TXRF
Berlin, Germany






TXRF – How does it work?

Principles of total reflection X-ray fluorescence (TXRF) spectroscopy



Total reflection X-ray fluorescence spectroscopy



Beam angle: $0^\circ / 90^\circ$

- Samples must be prepared on a reflective media
- Polished quartz glass or polyacrylic glass disc
- Dried to a thin layer, or as a thin film or microparticle

Principles of total reflection X-ray fluorescence spectroscopy




Samples for TXRF

- Powders: Direct preparation or as suspension
- Liquids: Direct preparation
- Always as a thin film, micro fragment or suspension of a powder
- Necessary sample amount: Low μg respectively μl range

Simple quantification

- ➔ Matrix effects are negligible due to thin layer
- ➔ Quantification is possible by internal standardization




Principles of total reflection X-ray fluorescence spectroscopy



In TXRF the samples are prepared as thin films or layers

- Matrix effects are negligible
- Quantification is possible




Principles of total reflection X-ray fluorescence spectroscopy




In TXRF the samples are prepared as thin films or layers

- Matrix effects are negligible
- Quantification is possible




- TXRF detects elements from Na(11) to U(92)
- The element sensitivities depend on the atomic number
- The sensitivity factors are calibrated ex works
- Quantification requires the addition of one standard element



Principles of total reflection X-ray fluorescence spectroscopy



- Mo tube offers lowest LLD
- W tube suitable for elements with $Z = 41$ to 53 (Nb to I)




The instrument S2 PICOFOX



Benchtop TXRF spectrometer S2 PICOFOX

- Metal-ceramic X-ray tube
 - Mo anode
 - Air-cooled
 - Optionally other tubes available
- Multilayer monochromator
- XFlash[®] silicon drift detector
 - Electro-thermally cooled
 - ≤ 149 eV @ MnK α 100 kcps
- Automatic version
 - 25 sample cassette





Wastewater analysis by TXRF spectrometry

Sample preparation

Suspensions and particulate matter



Suspensions can be analyzed right after dilution



Dilute sample
with distilled
water



Add internal
standard



Homogenize



Pipette on
carrier

Sample preparation


Suspensions and particulate matter



Dry through
heat/vacuum



Load the
instrument



Start data
acquisition

Industrial wastewater Cement




Situation

- Discharge of wastewater from a cement plant
- Threshold value of 500 $\mu\text{g/L}$ Pb to be fulfilled by law

Analytical task

- Filter water with fine particles
- Thickened sludge as a suspension




(cc) pd

Industrial wastewater Cement



Measurements and sample preparation



Industrial wastewater Cement



Results

- Filtered samples (with particles or clear) fulfill threshold values
- Raw sludge is close or above the limits filtration removes lead-containing solids
- The analysis of further elements is possible (not shown)

Pb ($\mu\text{g/L}$)	Filtered sewage	Filtered with particles	Sludge filtered	Sludge suspension
Conc. (sample 1)	n.d.	4.0	9.4	411
LLD (sample 1)	n.d.	2.5	2.0	9.5
Conc. (sample 2)	n.d.	n.d.	13	938
LLD (sample 2)	n.d.	n.d.	2.2	16

Industrial wastewater

Petro



Situation

- Control of wastewater in the petrochemical industry


Analytical task

- Analysis of wastewater
special interest in V, Mn, Ni, Hg



(cc) Leonard

Measurements and sample preparation



Results

- The elements V, Mn, Ni and Hg could be detected down low ppb
- The analysis of further elements is possible (not shown)

Conc. ($\mu\text{g/L}$)	Sample 1	Sample 2	Sample 3	Sample 4
V	< NWG	1.9	79	15
Mn	104	1.3	904	9.9
Ni	1.3	0.68	< NWG	5.4
Hg	< NWG	0.46	< NWG	4.0

Industrial wastewater

Oil/solvent mixtures




Situation

- Before discharging liquid mixtures of oils and solvents the content of heavy metals has to be verified

Analytical task

- All sewages have to be tested for the entire element range
- Therefore, all samples were analyzed by two instruments subsequently:
W unit for elements from Mo to Sb,
Mo unit for all other elements




(cc) Michael Meding

Industrial wastewater Oil/solvent mixtures



Measurements and sample preparation



Industrial wastewater

Oil/solvent mixtures



Results – sewage

- All elements from P to Hg were quantitatively analyzed
- Detection limits were < 10 ppb for most elements ($Z > 19$)
- Detection limits measured by the W unit were higher but fulfilled the required limits (e.g. Cd < 80 ppb, Sb < 160 ppb)

Results – oil samples and mixtures

- Best detection limits were in the range from 15 to 60 ppb (e.g. Br, Se, As, Zn, Ni, Cu)
- LLD measured by the W unit depend on sample type:
- Cd: 0.3 ppm in solvent, 1–2 ppm in oils and mixtures
- Sb: 0.6 ppm in solvent, 4 ppm in oils
- Mo: 0.2 ppm in solvent, 1 ppm in oils and mixtures

Improved sample preparation for achieving lower LLD is possible

Industrial wastewater TiO₂ production



Situation

- Before discharging the concentrations of several heavy metals they must be analyzed in order to fulfill legislative norms


Analytical task

- Wastewater has to be analyzed with regard to limit values for
 - Cr < 0,5 mg/l
 - Cd < 5,0 µg/l
 - Hg < 1,0 µg/l
- Concentrations in the mg/l range of additional elements (Fe, Ti and V) have to be supervised

Industrial wastewater TiO₂ production




Measurements and sample preparation



Industrial wastewater TiO₂ production




- Fast screening of wastewater delivers comparable results of TXRF and ICP-MS
- TXRF: Cd and Hg (conc. range 30 to 70 ppb) were not detectable due to line overlaps
- ICP-MS: S, Cl, Br, Rb, Sr were not detectable



Industrial wastewater TiO₂ production



- Detection limits



Water monitoring with bioindicators

Water monitoring with bioindicators

Introduction



Bioindicators are biological species used to ...

- monitor the health and integrity of an environment or ecosystem
- monitor for changes of metal pollutions in surface waters and sediments


General questions

- Accumulation through the nutrition chain?
- Impact of sublethal metal concentrations?
- Methods for fast screening?

Bioaccumulation of metals in the trophic chain



Element mass fraction increases at log scale



Metal screening of single Daphnia



Zooplankton Daphnia as bioindicator

- Describes biological availability of heavy metals, indicator for water quality
- Measurements of single organisms (< 100 μg) after selection by type, sex, growth status
- Trace element analysis of single organisms by ICP-OES and AAS impossible



Metal screening of single Daphnia



Sample preparation

- Rinsing of Daphnia
- Freeze drying
- Selection of individual organisms, weighing and size measurement
- Fix Daphnia on quartz carrier, add HNO_3 and internal standard
- Cold Plasma Ashing (2 h, 180°C)
no loss of volatile trace element compounds
- CPA reduces organic matrices and lowers background counts







TXRF measurement

- Typically 1000 s

Metal screening of single Daphnia



- Direct correlation of metal concentration and dry weight of bioindicator organism



Water monitoring with bioindicators


What are biofilms?



Dry sample
microscopy photo



Wet sample
underwater photo



Courtesy of UFZ Magdeburg

Significance of biofilms



Negative impact of biofilms

- Bioaccumulation of harmful chemicals in aquatic systems
- Biogenic corrosion of metals, concrete etc.
- Biofouling in pipeline systems
- Germination of water pipelines
- Contamination with pathogens in hospitals

Application of biofilms

- Self-purification of waters
- Binding and removal of toxic materials
- Biotechnological use
- Removal of organics and N-compounds in wastewater treatment plants


Biofilms

Sampling location



Elbe watershed

- Harz mountain range mining area since bronze age



Courtesy of UFZ Magdeburg



Biofilms

Sampling location



Sampling location "Schlüsselstollen"

- Gallery was built in 1879
- Drains the mining district at a length of 31 km



Courtesy of UFZ Magdeburg

Biofilms Sampling




Grown on round polycarbonate slide



Courtesy of UFZ Magdeburg

Biofilms

Sample preparation for TXRF



- Sampling and transport (4°C)
- Rinse with clean water
- Freeze drying
- Homogenization
- Weighing ca. 500 µg at Ultrabalance
- Digestion with $\text{HNO}_3/\text{H}_2\text{O}_2$
- Internal standard Sc
- 10 µl sample aliquot
- Drying on a hot plate 80°C

Analysis with TXRF

Accumulation factors in water and biofilms




Element	Accumulation factor water*	Accumulation factor biofilm*
Pb	1400	2200
Zn	630	18
Cu	11	140
Ni	8	7

*) relation of gallery to creek water

Biofilms: Pb > Cu > Zn > Ni > Fe > Ca

Sorption effects in biofilms



- Immediate accumulation of Pb and Ni during growth
- Release during longer exposition
- Slow accumulation of Cu and Zn

Biofilms

- Further investigations are planned on a river platform (meso cosmos) with “flow lines” simulating natural flow behaviours

Acknowledgement

- Margarete Mages, Helmholtz Centre for Environmental Research (UFZ) Magdeburg, Germany



Comparison with Atomic Spectroscopy

TXRF vs. AAS/ICP-OES



TXRF

- Small footprint/portable
- Microgram sample size
- %-ppb levels
- Single std. calibration
- Low maintenance
- Non-destructive
- Fast learning curve
- Fast sample prep

AAS/ICP-OES

- Large, fixed installation
- Milligram sample size
- ppm – sub ppb levels
- Multi-standard calibration
- High maintenance
- Destructive
- Long learning curve
- Laborious sample prep

Want ICP-MS? Bruker 820-MS



- Low DLs <1ppt for some elements in solution
- Wide range of elements (>75)
- Isotopic analysis capability



Innovative Ion optics 90° design



Ionen mirror


- Parabolic electrostatic field
- Reflects and focusses ions (*patented*)

Hollow structure of the ion mirror

- Photons and uncharged particles pass without contact

Arrangement of turbo pump

- optimal conditions for vacuum
- removes unwanted particles




Convenient technology

Collision Reaction Interface – CRI



Bruker's answer to common interferences?

- The Collision/Reaction Interface (CRI)
- H_2 gas injected to react with Ar^+ and ArX^+ molecular ions
- He gas used to separate other interferences
- Improved detection limits for several key elements (^{40}Ca , ^{56}Fe , ^{75}As , ^{80}Se)
- Increased productivity, most suitable for routine applications
- Fast and reproducible switching of CRI modes



Application labs in your neighborhood




Application labs with TXRF and ICP-MS at work or coming soon ...

- Berlin, Germany
- Coventry, UK
- Paris, France
- Milton, Canada
- Houston, Texas, USA
- Melbourne, Australia
- Seoul, Korea

Any Questions?

Please **type in** the questions you may have for our speakers in the **Questions form** and click **Submit**




Thank you for your attention!

Upcoming:

14th International Conference on
Total Reflection X-ray Fluorescence
and Related Methods
06–10 June 2011, Dortmund, Germany

www.bruker.com
www.s2picofox.com



www.bruker.com