

D4 ENDEAVOR

ANALYSIS OF SAMPLES CONTAINING FE WITH THE SOL-X DETECTOR

Introduction

Fluorescence radiation has an impairing effect on the quality of XRD measurement results if samples containing e.g. Fe or Mn are investigated using Cu-radiation. This is due to the significantly higher background and hence significantly lower peak-to-background ratio. Figure 1, for example, shows a measurement on a composite of Quartz and Hematite.

A diffracted beam monochromator is commonly used to suppress the fluorescence. The application of the monochromator drastically improves the peak-to-background ratio of the measurement (Fig. 2), but at the expense of a substantial loss of peak intensity. This eventually results in longer measurement times.

Switching to another X-ray wavelength is virtually the only alternative. Therefore, in cases where the loss of intensity can not be accepted from an analytical point of view, it is self-evident to have a diffracted beam set-up offering optimum peak intensities and a peak-to-background ratio.

This wish can now be fulfilled by the Sol-X energy dispersive detector – a new innovative XRD solution from Bruker AXS. It meets the requirements for a large X-ray active area as well as for an energy resolution to separate $\text{CuK}\beta$ radiation or $\text{FeK}\alpha$ fluorescence from $\text{CuK}\alpha$ simultaneously (Fig. 3). The user's main benefits are a reduction of the measurement time by a factor of about 3 and no extra effort when dealing with X-ray wavelength appropriate diffracted beam monochromators.

The Sol-X gives higher flexibility over diffracted beam monochromators, because it can be easily adjusted to the most commonly used wavelengths like $\text{CrK}\alpha$, $\text{CoK}\alpha$, $\text{MoK}\alpha$ and $\text{CuK}\beta$. With an active detector window area of $4 \times 15 \text{ mm}^2$ the entire beam of a long fine focus X-ray tube contributes to the XRD experiment, which results in higher count rates and, finally, smaller measurement times or higher quality of the data.

Instrumentation

- D4 ENDEAVOR
- 2.2 kW Cu long fine focus X-ray tube
- rotating sample stage
- 0.5° divergence slit system
- 0.5° anti-scatter slit system
- 4° Soller slits (primary and secondary)
- 0.1 mm receiving slit

Options:

- Ni-Filter, scintillation counter
- Graphite secondary monochromator, scintillation counter
- Sol-X, energy dispersive detector

Measurement parameters:

- Scan Mode: locked coupled
- Tube power: 40 kV, 35 mA
- Measurement range: 18° to 56° (2Theta)
- Step size: 0.02° (2Theta)
- Step time: 1 s




Figure 1: Measurement of a Hematite/Quartz mixture with Ni-Filter – high intensity, poor peak to background ratio. The green lines indicate the Quartz ICDD pattern, and the red lines indicate the Hematite pattern.




Figure 2: Measurement of a Hematite/Quartz mixture with secondary monochromator – low intensity, good peak to background ratio




Figure 3: Measurement of a Hematite/Quartz mixture with Sol-X solid state detector – high intensity, superior peak to background ratio

BRUKER AXS GMBH

BRUKER AXS, INC.

D-76187 KARLSRUHE
GERMANY

5465 EAST CHERYL PARKWAY
MADISON, WI 53711-5373
USA

TEL. (+49) (721) 595-2888
FAX (+49) (721) 595-4587
<http://www.bruker-axs.de>
Email: info@bruker-axs.de

TEL. (+1) (800) 234-XRAY
TEL. (+1) (608) 276-3000
FAX (+1) (608) 276-3006
<http://www.bruker-axs.com>
Email: info@bruker-axs.com