



with  
DAVINCI  
design

## The New D8 DISCOVER

- **Move up to the XRD<sup>2</sup> dimension in X-ray diffraction**
  - DAVINCI.SNAP-LOCK & DAVINCI.MODE: alignment-free component change with real-time recognition
  - TURBO X-RAY SOURCE: boost intensity for spot, line or micro focus applications
  - TWIST-TUBE: fast and easy switching from line to spot focus
  - VANTEC-500: add the extra dimension of XRD<sup>2</sup>
  - PATHFINDER optics: motorized switching between high-resolution and high-intensity beam paths

think forward

XRD