



D8 ADVANCE with DAVINCI.DESIGN

- **A series of innovations for the next era of X-ray diffraction**
 - TWIN/TWIN SETUP: Push-button switch between Bragg-Brentano and parallel-beam geometries
 - New DIFFRAC.SUITE with integrated virtual diffractometer
 - DAVINCI.SNAP-LOCK for alignment-free optics switch without tools
 - DAVINCI.MODE for real-time component recognition and configuration
 - TWIST-TUBE for fast and easy switching from line to point focus